

A.G. Scientific, Inc.



Announcing the AG Scientific Inhibitor Cocktail T-Shirt Contest!



Contest what Contest?

AG Scientific is excited about the interest and success researchers are having with AG's custom line of *Protease Inhibitor cocktails*. The use in proteomics is exploding and AG wants everyone to get on board! Researchers who are using tablet formed inhibitors will love AG's lyophilized inhibitors cocktail's instantaneous solubility, no chopping, no crunching, no waiting.

How to get into the contest:

Call or e-mail us and ask for a free sample of our Protease inhibitor cocktails and free t-shirt. We will ship it out to you. Then test drive AG's Protease Inhibitors, take a picture of you or your lab mates wearing the AG inhibitor T-shirt, and we will automatically enter into our raffle for an **Apple® - iPod shuffle® 2GB* MP3 Player (4th Generation - Latest Model) – Silver**. If you add a testimonial with your picture your chances double, we will put your name in twice. We will select a winner every 8 weeks so get a T-shirt and submit photos.

Kourosh Toumadje's Lab sporting a styling
AG Scientific Inhibitor Cocktail T-Shirt.
Submit your photo to: support@agscientific.com
Check the website for winners & updates
WWW.AGSCIENTIFIC.COM

AG Scientific Inc. ISO 9001:2008 Company
Phone 858-452-9925 Toll Free 877-452-9925
Email: Support@agscientific.com Web: www.agscientific.com

