

A.G. Scientific, Inc.

6450 Lusk Blvd. Suite E102

SAN DIEGO, CA

92121



CERTIFICATE of ANALYSIS

PRODUCT: **D-AP5**
[D(-)-2-Amino-5-phosphonovaleric acid]
C₅H₁₂N₀5P, M.W. 197.1, CAS# [79055-68-8]

PRODUCT NUMBER: A-1059

LOT NUMBER: B1308

APPEARANCE: White solid

SOLUBILITY: Soluble in water

ENANTIOMERIC PURITY: > 99%

PURITY: 99%

MELTING POINT: 233°C

APPLICATION NOTES: Potent NMDA antagonist. Active form of AP5. This product is not sterile.

STORAGE & HANDLING: Store at room temperature (+ 20°C) for up to 1 year.
Store solutions at -20°C for up to 3 month.
PROTECT FROM LIGHT!

CAUTION: For laboratory research & scientific manufacturing use only. Not for human or drug use. The pharmacological and toxicological properties of this product have not been fully investigated. Use caution when handling. Do not use this compound if you are not fully trained or are unaware of the hazards involved.

Verified: KT

A.G. Scientific, Inc, AN ISO 9001:2000 COMPANY
PHONE: (858) 452-9925, Toll FREE: (877) 452-9925, FAX: (858) 452-9926
Email: info@agscientific.com, Web site: www.agscientific.com