



CERTIFICATE OF ANALYSIS

PRODUCT:	Monoclonal Antibody to Bmf (9G10),
PRODUCT NUMBER:	B-1233
LOT NUMBER:	A9999
CROSSREACTIVITY SPECIES:	Human and Mouse
CLONE:	9G10
ISOTYPE:	Rat IgG2a
CONCENTRATION:	1 mg/ml
FORMULATION:	Liquid. Protein G-affinity purified antibody in PBS containing 0.02% sodium azide.
IMMUNOGEN:	Recombinant mouse Bmf.
SPECIFICITY:	Recognizes human and mouse Bmf.
APPLICATION:	Flow Cytometry Immunoprecipitation Western blot (excellent)

A.G. Scientific, Inc.

6450 Lusk Blvd. Suite E102

SAN DIEGO, CA
92121



STORAGE & HANDLING:

Store at -20°C. Avoid freeze/thaw cycles.

CAUTION: For laboratory research & scientific manufacturing use only. Not for human or drug use. The pharmacological and toxicological properties of this product have not been fully investigated. Use caution when handling. Do not use this compound if you are not fully trained or are unaware of the hazards involved.

Verified: KT

A.G. Scientific, Inc, AN ISO 9001:2000 COMPANY
PHONE: (858) 452-9925, Toll FREE: (877) 452-9925, FAX: (858) 452-9926
EMAIL: info@agscientific.com, WEB SITE: www.agscientific.com