



CERTIFICATE OF ANALYSIS

PRODUCT: **Monoclonal Antibody to Cystatin A
(human) (WR 23/2/3/3)**
[anti-Cystatin A Mab (human) (WR 23/2/3/3)]

PRODUCT NUMBER: C-1383

LOT NUMBER: A7986

CLONE: WR 23/2/3/3

ISOTYPE: Mouse IgG1 κ

IMMUNOGEN: Purified human cystatin A.

SPECIFICITY: Purified human cystatin A (Stefin A; ACPI).
The antibody was tested in several dilutions in a direct ELISA on plates coated with human cystatin A. Absorbance values fell below 50% of the maximal value at dilutions of 1 : 12,800. No cross-reactivity with cystatin B.

APPLICATIONS: In immunohistochemistry (paraffin sections); in Western blotting for the detection of human cystatin A.

PRESENTATION: Purified (Protein G)
200 μ g in 200 μ l PBS/0.01% Thimerosal

DILUTIONS: Working dilutions should be determined by titration assays.

STORAGE: Store at +4°C. Do not freeze!

CAUTION: For research use only. Not for human or drug use. The pharmacological and toxicological properties of this product have not been fully investigated. Use caution when handling. Do not use this compound if you are not fully trained or are unaware of the hazards involved.

Verified: 