



CERTIFICATE of ANALYSIS

PRODUCT: Polyclonal Antibody to NEDD8

PRODUCT NUMBER: S-1090

LOT NUMBER: A9414

SOURCE/HOST: From rabbit

IMMUNOGEN: Synthetic peptide corresponding to aa 20-32 (T²⁰DKVERIKERVEE³²) of human NEDD8.

SPECIFICITY: Recognizes human, mouse and rat NEDD8. Does Not cross-react with ubiquitin or SUMO.

FORMULATION: Liquid. Contains 50% glycerol containing 1mM sodium azide.

APPLICATIONS: Immunocytochemistry (1:500)
Immunohistochemistry
Immunoprecipitation
Western Blot (1:500)

STORAGE & HANDLING: Keep cool and dry at -20°C. Avoid freeze/thaw cycles.

CAUTION: For research use only. Not for human or drug use. The pharmacological and toxicological properties of this product have not been fully investigated. Use caution when handling. Do not use this compound if you are not fully trained or are unaware of the hazards involved.

Verified: KT