



PRODUCT DATA SHEET

PRODUCT: Polyclonal Antibody to Ubch9
[anti-Ubc9 (human) PAb]

LOT NUMBER: A7400

PRODUCT NUMBER: S-1092

HOST: From Sheep.

CONCENTRATION: 0.5 mg/ml

SPECIFICITY: Recognizes Ubch9.

SPECIES

CROSSREACTIVITY: Human

IMMUNOGEN: Synthetic peptide corresponding to aa 147
158 (R¹⁴⁷VRAQAKKFAPS¹⁵⁸) of C-terminal Ubch9.

FORMULATION: Liquid. Affinity purified antibody in 100 mM TRIS-HCl, pH 7.5,
containing 10% glycerol and 0.01% sodium azide.

APPLICATION: Immunocytochemistry, Western Blot (1:3'000)

STORAGE & HANDLING: Keep cool and dry at -80°C. Once thawed add glycerol and store
in a frost free freezer at -20°C for up to 6 months. Avoid
freeze/thaw cycles.

CAUTION: For research use only. Not for human or drug use. The pharmacological and toxicological properties of this product have not been fully investigated. Use caution when handling. Do not use this compound if you are not fully trained or are unaware of the hazards involved.

Verified: 